



# EMS Physio Ltd

## Variable Height Tilt Table

The Standard Tilt Table has been designed for a wide range of rehabilitation procedures. Featuring electric height and tilt operation to a full 90° for patient “walk off” procedures. The table has a safe working load of 250Kg and is supplied with adjustable height removable handgrips for added security. Angle indicators are fitted on both sides for reference and large braked wheels enable the table to be moved and positioned with ease. A wide choice of options are available and include dual adjustable foot boards, divided leg design with pommel and a “see-through” Perspex work table with angular adjustment.




### Standard Features

- UK designed & manufactured
- Lifetime frame & 5 year parts guarantee
- Hand switch control of height and tilt operations
- Fixed single footboard as standard
- Pair of height adjustable hand grips
- Angle inclinometers on both sides
- Comfort foam upholstery
- Hard wearing antimicrobial vinyl
- Tilt facility to 90°
- Set of three patient harnesses
- Large, independently braked wheels
- Wide choice of options available

### Overall Dimensions

Length: 200cm (78”) approx  
 Width: 85cm (33”) approx  
 Height: 55cm - 97cm approx

### Upholstery Dimensions



### Standard Features



Foot Board



Angle Indicators



Large Braked  
Wheels



Patient Hand Grips




Patient Harnesses

NHS Approved Supplier



t: +44 (0) 1235 772272  
 e: sales@emsphysio.co.uk  
 www.emsphysio.co.uk

## Operational Parameters



EMS Physio Ltd order codes  
COU007 Tilt Table

## For Your Piece of Mind

- Lifetime frame and 5 year parts (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average life span of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2015 by LQRA

## Standards

- IP Rating - IPX6
- EE Classification - EU MDR 2017/745 Annex VIII Class 1
- Electrical Standard - IEC 60601-1
- Vinyl flame retardant to CRIB 5 (BS 5852)
- Foam flame retardant to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

## Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required.

## Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note; Unless specified Sapphire upholstery will be supplied

## Optional Accessories



Divided Leg with Pommel



Divided Foot Board



Spare Harness Set



Perspex Patient Table



Head Cushion



Battery Backup



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond

**EMS PHYSIO Ltd**  
Downsview Road  
Wantage  
Oxfordshire OX12 9FE

**t: +44 (0) 1235 772272**  
**f: +44 (0) 1235 763518**  
**e: sales@emsphysio.co.uk**  
**www.emsphysio.co.uk**