



A very strong and rigid 2 section couch designed for various rehabilitation procedures, general examinations and massage etc. This model requires minimal maintenance and is suitable for constant heavy usage. The couch is available with a choice of hydraulic or electric height adjustment and features a manually adjustable backrest for supporting the patient in a sitting position, or it can be positioned as required for prone or supine lying procedures.

A wide range of options are available to personalise the couch to individual user requirements.



Standard Features

- UK manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Safe working load & lift capacity 250kg
- Electric or hydraulic operation
- Low running costs
- Anti-entrapment gas strut
- Hardwearing antimicrobial vinyl
- Retractable wheel base
- Adjustable foot for uneven floors
- Comfort foam
- Easy operation
- Hoist clearance 13cm (5")

Overall Dimensions

Length: 188cm (74") approx
 Width: 63cm (25") approx
 Height: 47cm - 97cm approx

Upholstery Dimensions




Options & Accessories



Breathing Hole



Paper Roll Holder



Side Security Rails



D Cushion



Footswitch



Motorised Backrest

Operational Parameters



Couch with backrest negative and optional paper roll holder



Couch at minimum height with optional paper roll holder

EMS Physio Ltd order codes
COU001 Hydraulic with breathing slot
COU002 Electric with breathing slot

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note; Unless specified Sapphire upholstery will be supplied



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond

For Your Piece of Mind

- Lifetime frame and 5 year parts (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average life span of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2015 by LQRA

Standards

- IP Rating - IPX6
- EE Classification - EU MDR 2017/745 Annex VIII Class 1
- Electrical Standard - IEC 60601-1
- Vinyl flame retardant to CRIB 5 (BS 5852)
- Foam flame retardant to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required.

EMS PHYSIO Ltd
Downsview Road
Wantage
Oxfordshire OX12 9FE

t: +44 (0) 1235 772272
f: +44 (0) 1235 763518
e: sales@emsphysio.co.uk
www.emsphysio.co.uk