



EMS PHYSIO LTD

Manufacture

Part Number: 2702 Couch 2 section

User Manual/Specification

Portable Couch with carry case

Weighing in at only 14.5kg this aluminium couch represents the perfect solution for the mobile massage therapist while retaining great strength from the frames extensive use of aluminium. Free carry bag included.

Arm and face cradle optional extra – only available in blue or white upholstery



Description/Features

- Extremely light & strong aluminium frame
- Multi Adjustable Table Height – push button height control
- 10 position flip up back rest section
- Safety Lock Support System
- Steel Cable Tension System
- Set Up Protectors
- High Density Dual Layer Foam 5cm
- Rounded Corners That Are Less Prone To Damage
- Carry Bag with Shoulder Strap & Pockets



Specification

Extended therapy table size (exclude head and arm rest)

Length: 182 cm

Width: 60 cm

Height: 66 cm – 86 cm


Vinyl: High quality PU, oil and waterproof

Vinyl Thickness: 4cm thick made up of dual layers, one soft layer and one thicker layer total depth 5cm
Colour: Black, white, Navy blue: Aluminium frame

Working load: 204kg

Table weight: (Table Only): 14.5 kg

Size when folded: Length: 91.4cm x Depth:15.2cm x Height: 60.96cm



Ordering information:

Part number: 2702

sales@emsphysio.co.uk

+44 (0) 1235 772272

www.emsphysio.co.uk

EMS Physio Ltd, Grove Business Park, Downsview Road, Wantage, Oxfordshire,
OX12 9FE

Company Registration Number: 601612